

# Passive sampling of nonpolar compounds in sediments

Kees Booij (NIOZ, NL) Foppe Smedes (Deltares, NL)

SERDP/ESTCP Passive sampling session 2010-12-02

#### **Report Documentation Page**

Form Approved OMB No. 0704-0188

Public reporting burden for the collection of information is estimated to average 1 hour per response, including the time for reviewing instructions, searching existing data sources, gathering and maintaining the data needed, and completing and reviewing the collection of information. Send comments regarding this burden estimate or any other aspect of this collection of information, including suggestions for reducing this burden, to Washington Headquarters Services, Directorate for Information Operations and Reports, 1215 Jefferson Davis Highway, Suite 1204, Arlington VA 22202-4302. Respondents should be aware that notwithstanding any other provision of law, no person shall be subject to a penalty for failing to comply with a collection of information if it does not display a currently valid OMB control number.

1. REPORT DATE 02 DEC 2010	2. REPORT TYPE	3. DATES COVERED <b>00-00-2010 to 00-00-2010</b>	
4. TITLE AND SUBTITLE	5a. CONTRACT NUMBER		
Passive sampling of nonpolar compounds in sediments		5b. GRANT NUMBER	
	5c. PROGRAM ELEMENT NUMBER		
6. AUTHOR(S)		5d. PROJECT NUMBER	
	5e. TASK NUMBER		
	5f. WORK UNIT NUMBER		
7. PERFORMING ORGANIZATION NAME(S) AND AI Royal Netherlands Institute for Sea RoNetherlands,	8. PERFORMING ORGANIZATION REPORT NUMBER		
9. SPONSORING/MONITORING AGENCY NAME(S) AND ADDRESS(ES)		10. SPONSOR/MONITOR'S ACRONYM(S	
		11. SPONSOR/MONITOR'S REPORT NUMBER(S)	

12. DISTRIBUTION/AVAILABILITY STATEMENT

Approved for public release; distribution unlimited

13. SUPPLEMENTARY NOTES

Presented at the 15th Annual Partners in Environmental Technology Technical Symposium & Workshop, 30 Nov ? 2 Dec 2010, Washington, DC. Sponsored by SERDP and ESTCP. U.S. Government or Federal Rights License

14. ABSTRACT

Passive sampling is a powerful tool for determining the concentrations of truly dissolved contaminants in pore waters as well as for measuring the size of the exchangeable contaminant fraction in the sediment phase. Application of passive sampling of nonpolar contaminants in sediments is technically easy, but conceptually difficult, compared with traditional batch water sampling. The sorption capacity of the samplers and the (equivalent) water sampling rate are the two parameters that are needed to conceptualize the sampling process. These parameters also determine the value of the detection limits and the rate at which passive samplers equilibrate with their environment. After discussing contaminant uptake by passive samplers in an infinite water volume, their behavior in sediments will be dealt with using two practical examples. Exposure to stagnant sediments will be contrasted with exposure to sediment slurries. Passive sampler exposure to sediment slurries at low sampler/sediment phase ratios allows concentrations of dissolved contaminants in the pore water to be estimated. Exposure at high sampler/sediment phase ratios yields information on the size of the contaminant fraction in the sediment that is easily exchangeable with the pore water. Guidelines for designing passive sampler exposures to contaminated sediments will be discussed, and sources of uncertainties in the results will be highlighted.

15. SUBJECT TERMS								
16. SECURITY CLASSIFICATION OF:			17. LIMITATION OF ABSTRACT	18. NUMBER OF PAGES	19a. NAME OF RESPONSIBLE PERSON			
a. REPORT unclassified	b. ABSTRACT unclassified	c. THIS PAGE unclassified	Same as Report (SAR)	30	1.00.01.00.00.00.00.00.00.00.00.00.00.00			

#### PASSIVE SAMPLING OF NONPOLAR COMPOUNDS IN SEDIMENTS

KEES BOOIJ
Royal Netherlands Institute for Sea Research
P.O. Box 59
Texel, 1790 AB NETHERLANDS
+31 222 369 463
kees.booij@nioz.nl

CO-PERFORMER: Foppe Smedes (Deltares)

Dassive sampling is a powerful tool for determining the concentrations of truly dissolved **I** contaminants in pore waters as well as for measuring the size of the exchangeable contaminant fraction in the sediment phase. Application of passive sampling of nonpolar contaminants in sediments is technically easy, but conceptually difficult, compared with traditional batch water sampling. The sorption capacity of the samplers and the (equivalent) water sampling rate are the two parameters that are needed to conceptualize the sampling process. These parameters also determine the value of the detection limits and the rate at which passive samplers equilibrate with their environment. After discussing contaminant uptake by passive samplers in an infinite water volume, their behavior in sediments will be dealt with, using two practical examples. Exposure to stagnant sediments will be contrasted with exposure to sediment slurries. Passive sampler exposure to sediment slurries at low sampler/sediment phase ratios allows concentrations of dissolved contaminants in the pore water to be estimated. Exposure at high sampler/sediment phase ratios yields information on the size of the contaminant fraction in the sediment that is easily exchangeable with the pore water. Guidelines for designing passive sampler exposures to contaminated sediments will be discussed, and sources of uncertainties in the results will be highlighted.



# Nonpolar samplers for analytes with $log K_{ow} > 4$

Technically simple
Conceptually difficult (initially)





#### Theory

- working principles
- concepts & link to batch sampling
  - sampling rates and sorption capacity
  - o equilibration times and detection limits
- sampling rate calibration
- finite water volumes and sediment slurries


#### Practice

- pore water extraction (measuring  $C_{\rm w}$ )
- sediment extraction (measuring accessible fractions)
- Achievements & caveats




# PSDs in water (infinite volume)

## Workingprinciple



## Uptake kinetics



## Linear uptake stage

$$N = R_{\rm S} C_{\rm W} t$$

$$\frac{N}{C_w}$$
 = effectively extracted water volume =  $R_s t$ 

$$N \neq R_s t C_w$$

amount = water volume x concentration

time

amount



#### Equilibrium stage

$$\frac{N}{m_{\rm s}} = K_{\rm sw} C_{\rm w}$$


$$\frac{N}{C_{\rm w}}$$
 = effectively extracted water volume =  $m_{\rm s} K_{\rm sw}$ 

$$N \neq (m_s K_{sw}) C_{w}$$

amount = water volume x concentration




## Uptake model (infinite water volume)





### Equilibration time (TWA window)





#### **Detection limits**


Linear uptake stage

$$V_{\text{wLUS}} = R_{\text{s}} t$$

take the lowest of these

Equilibrium stage


$$V_{\text{wEqS}} = m_{\text{s}} K_{\text{sw}}$$



$$DL_{CW} = \frac{DL_N}{V_W}$$



#### Uptake rate control



membrane-control

$$R_{\rm s} \sim K_{\rm sw} D_{\rm m}$$

**WBL-control** 


$$R_{\rm s} \sim D_{\rm w}^{2/3}$$



#### Sampling rates $(R_s)$

compound properties (diffusion coefficient, hydrophobicity)

- temperature
- flow rate
- suspended solids concentration






#### Calibration of uptake kinetics

#### Performance reference compounds

- not occurring in the water
- spiked into PSD before exposure
- dissipation follows:  $\frac{N_s}{N_0} = \exp\left(-\frac{R_s}{K_{sw}} \frac{t}{V_s}\right)$






# PSDs in water (finite volume)




#### Infinite water volume



"infinite" means:  $V_{\rm w} >> m_{\rm s} K_{\rm sw}$ 



#### Finite water volume



"finite" means:  $V_{\rm w} \cong m_{\rm s} K_{\rm sw}$ 


phenanthrene in 10 L water 1 g silicone  $m_s K_{sw} = 13 L$ 

#### faster equilibrium

- at a lower level
- because of depletion



### Sampling in sediment slurries



Pore water sniffing:  $m_{\text{sed}} K_{\text{d}} + V_{\text{w}} >> m_{\text{s}} K_{\text{sw}}$  $C_{\text{w}}(t) \cong C_{\text{w,0}}$ 

## Example 1: pore water sniffing

PAHs/PCBs/HCB in harbor sediments

ES&T 37: 4213-4220,2003



### LDPE strips in sediment slurries and stagnant sediments

0

10

20


30

*t* (d)

40

50


slurry



 $m_{\rm s}$   $K_{\rm sw} \approx 300 \, \text{L} \, (\oplus \text{ or } \odot)$ 







8 - 7 - 2° y 60 6 - 5 - 6 7 8 log K<sub>ow</sub>

from  $C_{\text{sed}}/C_{\text{w}}$ 

from uptake rates in stagnant sediments

$$N \sim A C_{\rm w} \sqrt{K_{\rm d}} \sqrt{t}$$

## Example 2: nondepletive and depletive extraction

PAHs/PCBs/HCB in estuarine sediments

Smedes, unpublished results

#### Scaling the experiment

- ~20 g sediment (2.5% oc)
- ~200 mL water

0.2, 2, 20 g silicone sheet


#### 48d exposure

$$m_{\rm sed} K_{\rm d}$$
  $\approx$  120 L  $m_{\rm s} K_{\rm sw}$   $\approx$  10, 100, 1000 L (©) sniffing depleting



#### High sampler mass ⇒ depletion


### 20 g dw sediment, $f_{oc}$ =0.025






# Equilibrate samplers more mass $\Rightarrow$ lower $C_w \Rightarrow$ lower $C_{sed}$

#### Fluoranthene





Indeno[1,2,3-cd]pyrene



accessible : 14%

inaccessible: 86%

# PSDs for sediment quality assessment: achievements & caveats

#### **Achievements**

- use  $R_s$  and  $m_s K_{sw}$  to link PS to batch sampling
- non-depletive  $\rightarrow C_w$  (need  $R_s$  calibration)
- depletive → accessible fraction (demonstrate equilibrium)
- scaling rules are straightforward ( $m_{\text{sed}} \& m_{\text{sampler}}$ )

#### **Caveats**

- availability of reliable K<sub>sw</sub> values
- among-manufacturer variability of polymer properties
- limited PSD intercomparison





